

ULTRASONIC ECHOSCOPE GS200/GS200i

Reflection and Transmission Measurements

A-scan, B-scan, TM-mode, CT-mode & B-imaging



User Manual



Revision: Version 0.70 EN, Januar 26, 2016

Gesellschaft für Angewandte Medizinische Physik und Technik mbH
Hallesche Straße 99F
06217 Merseburg
Phone: +49 3461/278691-0
Fax: +49 3461/278691-101
E-mail: info@gampt.de

Visit us on the internet: www.gampt.de

© GAMPT mbH

All rights are reserved for GAMPT mbH. No part of this manual may be reproduced or processed, copied or distributed in any form without the prior approval of GAMPT mbH.

GAMPT mbH is not liable for damages occurring as a result of incorrect use, or for repairs and alterations carried out by third parties not authorized by GAMPT mbH.

Subject to technical alterations. Errors excepted.

Printed: 02/01/2016

Version: 0.70 EN

Table of Contents

1	Safety information	4
2	Scope of delivery	4
3	Introduction	5
4	Ultrasonic echoscope	6
4.1	The basic unit GS200	6
4.2	The extended device version GS200i	6
4.3	Control elements on the front side	7
4.3.1	Receiver/transmitter unit (RECEIVER/TRANSMITTER)	8
4.3.2	TGC unit (Time Gain Control)	8
4.3.3	BNC outputs (MONITOR)	11
4.4	Connection ports on the rear side of GS200/GS200i	12
4.4.1	Fuse replacement	12
4.4.2	External triggering	13
4.5	Ultrasonic probes	14
4.5.1	Single-element transducers (GS200/GS200i)	14
4.5.2	Array transducer (GS200i)	15
5	Software	17
5.1	General program structure	17
5.1.1	User interface of the program	17
5.1.2	Function switches	18
5.1.3	Diagrams	18
5.1.3.1	Diagram tool bars	18
5.1.3.2	Zooming and shifting of diagram contents using the mouse	19
5.1.3.3	Measuring cursors	22
5.1.4	Print command	23
5.2	Program start	24
5.2.1	Initialization.....	24
5.2.2	Simulation mode.....	25
5.3	Echoscope operation with single-element transducers (GS200/GS200i)	26
5.3.1	Starting, interrupting, and completing the measuring mode	26
5.3.2	Visualization of the measurement signals in an A-scan diagram.....	26

5.3.3	Measuring and program parameters	27
5.3.3.1	Params tab	28
5.3.3.2	USB tab	31
5.3.3.3	Draw tab	32
5.3.3.4	Scanner tab	35
5.3.4	Information on device settings and program status.....	36
5.3.4.1	Device information (Front panel)	36
5.3.4.2	Information on program and measuring operation	37
5.3.4.3	Scanner status	40
5.3.5	A-mode.....	41
5.3.5.1	Time-of-flight and amplitude measurements	42
5.3.5.2	Measuring accuracy	44
5.3.5.3	Frequency spectrum and frequency filter	45
5.3.5.4	Cepstrum and cepstrum filter	48
5.3.6	B-mode.....	51
5.3.7	M-mode.....	54
5.3.8	CT-mode	56
5.3.8.1	Ultrasonic computer tomography	56
5.3.8.2	CT-mode tab.....	57
5.3.8.3	Internal and external triggering of CT-scans.....	59
5.3.8.4	Procedure of a CT measurement	61
5.4	Echoscope operation with array transducer (GS200i).....	64
5.4.1	Characteristic features of the Imaging Mode.....	64
5.4.2	Start, interrupt and terminate the measuring mode	65
5.4.3	Scan modes	65
5.4.4	Measuring and program parameters	66
5.4.4.1	Main tab	67
5.4.4.2	TGC tab.....	69
5.4.4.3	M-Line tab	70
5.4.4.4	Image tab.....	70
5.4.4.5	Palette tab	73
5.4.4.6	Cine tab.....	74
5.4.5	Features of the info field.....	75
5.4.5.1	Probe detection	75
5.4.5.2	Measuring cursor	76
6	Accessories	77






6.1	CT equipment	77
6.1.1	CT control unit UCT200.....	77
6.1.2	CT scanner.....	78
7	Experiments	80
8	Technical data	81
8.1	Echoscope.....	81
8.2	Ultrasonic probes	82
8.3	CT devices.....	82
9	Literature	83
10	Revisions of the user manual	83
11	Index	84



Please read through the whole user manual before using the device and accessories.

1 Safety information

Prior to taking the ultrasonic echoscope and the accessories into operation, please carefully read the notes below to ensure your own health and safety and operational safety of the device.

-  *The opening slits on the device are for ventilation and must be kept absolutely free in order to prevent the device from overheating. It is recommended to use the folding feet on the device.*
-  *Make sure that the voltage values and fusing stated for the device are adhered to in the electricity supply.*
-  *Never try to insert objects through the openings on the device, because this can lead to short circuits or electric shocks.*
-  *Only the ultrasonic transducers supplied by GAMPT mbH must be connected to the ports named "PROBE 1" and "PROBE 2". Caution! During transmission mode voltage up to 300 V may be applied to the ports.*
-  *Attention, the ultrasonic echoscope is a laboratory device not a medical product! The ultrasonic probes must not be used with persons or other living beings.*



Compliance with the special operating and safety information in the individual chapters is mandatory.

2 Scope of delivery

Device:	GS200	GS200i
Order number:	10400	10410
Connecting cable:	<ul style="list-style-type: none"> • Mains cable • USB cable 	<ul style="list-style-type: none"> • Mains cable • USB cable
Measurement software:	<ul style="list-style-type: none"> • GS-EchoView 	<ul style="list-style-type: none"> • GS-EchoView
Ultrasonic transducer:		<ul style="list-style-type: none"> • Convex array transducer

